





## SAFETY DATA SHEET

### SECTION 01 - PRODUCT AND COMPANY IDENTIFICATION

Product Identifier ANISOLE		Product Use Laboratory use	
Chemical formula C <sub>6</sub> H <sub>5</sub> OCH <sub>3</sub>		Product code AR-0929	Molar weight 108,14
Chemical name / Commercial name / Synonymous ANISOLE ANHYDROUS, METHOXYBENZENE, METHYL PHENYL ETHER, PHENYL METHYL ETHER			
Supplier's name Laboratoire MAT		Address-Street 610, Adanac Street	
City Québec		Province Québec	
Postal code G1C 7B7	Internet www.labmat.com	Phone number 418-660-8666 Lun-Ven 8h-16h	
Emergency phone	418-660-8666 Lun-Ven 8h-16h		CENTRE ANTI-POISON DU QUÉBEC 800-463-5060
Date SDS 3/18/2024	SDS Prepared by Laboratoire MAT		E-Mail labmat@labmat.com

## SECTION 02 - HAZARDS IDENTIFICATION

<b>Classification WHIMS / GHS</b>	Skin corrosion/irritation - Skin irritation category 2 Flammable liquids category 3 Specific target organ toxicity - Single exposure category 3 Serious eye damage/eye irritation - Eye irritation category 2
<b>Signal Word</b>	WARNING
<b>Hazards statements (H)</b>	H226 Flammable liquid and vapor. H315 Causes skin irritation. H319 Causes serious eye irritation. H336 May cause drowsiness or dizziness.
<b>Precautionary statements (P)</b>	P210 Keep away from heat/sparks/open flames/hot surfaces. — No smoking. P233 Keep container tightly closed. P240 Ground/bond container and receiving equipment. P241 Use explosion-proof electrical/ventilating/lighting equipment. P242 Use only non-sparking tools. P243 Take precautionary measures against static discharge. P261 Avoid breathing dust / fume / gas / mist / vapors / spray. P264 Wash the areas of the body that have been in contact with the product after handling. P271 Use only outdoors or in a well-ventilated area. P280 Wear protective gloves/protective clothing/eye protection/face protection. P302 + P352 IF ON SKIN: Wash with plenty of soap and water. P303 + P361 + P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. P304 + P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P312 Call a POISON CENTER or doctor/physician if you feel unwell. P321 Specific treatment (see section 4 of the SDS and on this label). P332 + P313 If skin irritation occurs: Get medical advice/attention. P337 + P313 If eye irritation persists: Get medical advice/attention. P362 + P364 Take off contaminated clothing and wash it before reuse. P370 + P378 In case of fire: Use water spray or alcohol-resistant foam, or dry powder or carbon dioxide for extinction. P403 + P233 Store in a well-ventilated place. Keep container tightly closed. P403 + P235 Store in a well-ventilated place. Keep cool. P405 Store locked up. P501 Dispose of contents and container in accordance with local, regional and national regulations, or contact a specialist waste disposal company.
<b>PICTOGRAMS</b>	
<b>Other dangers</b>	NFPA (Risk: 0=No risk; 1=Slight; 2=Moderate; 3=Signifiant; 4=Extreme)
	<b>Health</b> 2 <b>Fire</b> 2 <b>Reactivity</b> 0 <b>Special danger</b>

## SECTION 03 - COMPOSITION/INFORMATION ON INGREDIENTS

Ingrédients (Dénomination chimique / synonymes)	Numéro CAS et tout identificateur unique	Concentration (%P/P)
Anisole	100-66-3	<=100

## SECTION 04 - FIRST AID MEASURES

<b>Eye contact</b>	Wash eyes with large amounts of water for at least 15 minutes while holding eyelids apart to rinse eyes. If irritation persists, seek medical attention.
<b>Skin contact</b>	Wash skin with plenty of water for at least 15 minutes. Remove soiled clothing. If irritation persists, seek medical attention.
<b>Inhalation</b>	If breathed in, move person into fresh air. If breathing is difficult, give oxygen. Consult a physician.
<b>Ingestion</b>	If the person is conscious, give water to drink. Never give anything by mouth to an unconscious person. Consult a physician.
<b>Most important symptoms and effects (acute and delayed)</b>	Main symptoms of high exposure: Skin lesions. Headaches. Vertigo. Nausea and vomiting. Tiredness. Ref. section 11.
<b>Immediate medical attention and special treatment, if necessary</b>	In case of medical consultation, keep this sheet available.
<b>General advice</b>	Show this safety data sheet to the doctor in attendance.

## SECTION 05 - FIREFIGHTING MEASURES

<b>Flammability</b>	Yes
<b>Suitable extinguishing media</b>	Water spray, carbon dioxide, dry chemical powder, alcohol or polymer foam.
<b>Unsuitable extinguishing media</b>	Data not available.
<b>Hazardous combustion products</b>	Hazardous combustion products formed under fire conditions: Carbon oxides.
<b>Specific hazards of the dangerous product</b>	Steam can travel a great distance and ignite on sources of ignition such as heaters, electrical appliances, cigarettes, sparks, etc. Containers exposed to fire may explode. May react violently with incompatible products (Ref Section 10).
<b>Special protective equipment and precautions for firefighters</b>	Discard incompatible substances if this can be done without risk. Firefighters should be equipped with standard protective equipment, fireproof clothing, face mask, gloves, protective boots and, where appropriate, self-contained breathing apparatus.

## SECTION 06 - ACCIDENTAL RELEASE MEASURES

<b>Methods and materials for containment and cleaning up / Personnel precautions, protective equipment</b>	Evacuate personnel to safe areas. Cut off all sources of ignition. Avoid the accumulation of charges electrostatic. Avoid breathing vapors, mist or gas. Absorb the product with sand or vermiculite. Dilute residues with water, clean and rinse. Ensure a good ventilation of the premises. Dispose of residues in a container for disposal of hazardous materials. When handling, wear suitable safety equipment. Use breathing apparatus if necessary. Vapors can accumulate in low areas. Do not let product enter drains.
--	---

## SECTION 07 - HANDLING AND STORAGE

<b>Conditions for safe storage</b>	Store in a cool, dry place. Keep container tightly closed and store away from heat, moisture, combustible and incompatible products. Keep away from sources of ignition - No smoking. Take measures to prevent the accumulation of electrostatic charges. Use venting and electrical equipment that is grounded and does not produce ignition sources (sparks). Protect from the sun's rays. Protect from the sun's rays.
<b>Methods of handling</b>	Always open containers slowly to allow any excess pressure to vent. Avoid ingestion and inhalation. Ensure good ventilation. Wear personal protective equipment when handling. Always ensure good ventilation. Transport according to TDG (ref Section 14)

## SECTION 08 - EXPOSURE CONTROLS/PERSONAL PROTECTION

### Workplace control parameters

Components	CAS-No.	Control parameters	Value	Basis
Anisole	100-66-3	TLV, TWA, STEL	No data available	Canada. Alberta, Occupational Health and Safety Code (table 2: OEL)
		TLV, TWA, STEL	No data available	Canada. British Columbia OEL
		TLV, TWA, STEL	No data available	Québec. Regulation respecting occupational health and safety, Schedule 1, Part 1: Permissible exposure values for airborne contaminants

<b>Ventilation</b>	Fan.
<b>Respiratory</b>	If the permissible levels are exceeded, use a mechanical filter / cartridge against NIOSH vapors or a respirator with air supply.
<b>Gloves</b>	Handle with gloves.
<b>Eyes</b>	Safety goggles with safety shutters.
<b>Shoes</b>	Safety shoes.
<b>Clothing</b>	Labcoat. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.
<b>Engineering control</b>	Have safety showers and eyewash stations in the workplace in case of an emergency and a ventilation system to maintain the level of concentrations in the air below the exposure limit values.

## SECTION 09 - PHYSICAL AND CHEMICAL PROPERTIES

<b>Physical state</b>	Liquid.
<b>Appearance</b>	Incolore à jaune.
<b>Odour</b>	Douce et aromatique, odeur anisée.
<b>Odour threshold</b>	Data not available
<b>pH</b>	Sans objet.
<b>Melting point / Freezing point</b>	-37°C
<b>Initial boiling point</b>	154°C
<b>Boiling range</b>	Data not available
<b>Flash point</b>	51°C
<b>Evaporation rate</b>	Data not available
<b>Flammability</b>	Yes
<b>Lower flammable / Explosive limit</b>	0.34%
<b>Upper flammable / Explosive limit</b>	6.3%
<b>Vapour pressure</b>	10.0 mm @ 42.2°CmmHg
<b>Vapour density</b>	3.72-
<b>Relative density</b>	0.995g/cm <sup>3</sup>
<b>Solubility</b>	Insoluble dans l'eau. Miscible dans l'alcool et l'éther.
<b>Partition coefficient water/n-octanol</b>	Data not available
<b>Auto-ignition temperature</b>	475°C
<b>Decomposition temperature</b>	Data not available
<b>Viscosity</b>	1.52 cP @ 15°C cps.

## SECTION 10 - STABILITY AND REACTIVITY

<b>Reactivity</b>	May react violently with incompatible substances. Flammable product, may ignite with source of ignition, if temperature above flash point.
<b>Chemical stability</b>	Stable under recommended storage conditions. Sensitive to heat.
<b>Possibility of hazardous reactions</b>	In case of strong heating, possibility of formation of explosive mixtures with air. May react violently with incompatible substances. Flammable liquid and vapors. May form explosive peroxides.
<b>Conditions to avoid, including static discharge, shock or vibration</b>	This product quickly absorbs moisture from the air. Avoid contact with incompatible materials and extreme temperatures. Reacts with air to form peroxides. Avoid the accumulation of static electricity.
<b>Incompatible materials</b>	Strong oxidizing agents (nitric acid, perchloric acid, peroxides, chlorates and perchlorates), strong acids, heat and moisture.
<b>Hazardous decomposition products</b>	Toxic vapors of carbon monoxide and dioxide.

## SECTION 11 - TOXICOLOGICAL INFORMATION

### ANISOLE

<b>Routes of exposure</b>	Ingestion and skin contact.
<b>Acute exposition effects / symptoms:</b>	By exposure route below.
<b>- Eyes</b>	Irritation and tearing.
<b>- Skin</b>	Irritation and dermatitis. It can cause redness, flaking and cracking.
<b>- Inhalation</b>	Irritation of the mucous membranes and respiratory tract. Nervous disorders, chest pain, cough, dyspnea, headache, dizziness, drowsiness, nausea and vomiting.
<b>Acute toxicity (Ingestion)</b>	Irritation of the mucous membranes. Gastrointestinal disorders, liver damage, cramps, diarrhea, headache, dizziness, drowsiness, seizures, nausea and vomiting.
<b>Chronic exposure effects / symptoms</b>	Burning sensation, dermatitis, nervous disorders, chest pain, cough, dyspnea, laryngitis, headache, dizziness, drowsiness, irritability, muscle weakness, weight loss and loss of appetite, nausea and vomiting.
<b>DL50 (specify species and route of entry)</b>	LD50 Oral - Rat - 3700 mg/kg. LD50 Dermal: Data not available.
<b>CL50 (specify species and route of entry)</b>	CL50 inhalation - Rat 6.5 mg/L - 4 h

## SECTION 12 - ECOLOGICAL INFORMATION

<b>Ecotoxicity</b>	Toxicity to fish: LC50 - Danio rerio (zebra fish) - > 1 mg/l - 96 h. Toxicity to daphnia and other aquatic invertebrates: Daphnia magna (Water flea) - 11.05 mg/l - 24 h.
<b>Persistence and degradability</b>	Biodegradability aerobic - Exposure duration 14 days: 56%. Readily biodegradable.
<b>Bioaccumulative potential</b>	Data not available.
<b>Mobility in soil</b>	Data not available.
<b>Other adverse effects</b>	Data not available.

## SECTION 13 - DISPOSAL CONSIDERATIONS

<b>Waste Disposal Method</b>	Dispose of contents and container in accordance with local, regional and national regulations, or contact a specialist waste disposal company.
<b>Contaminated Packaging</b>	Dispose of as unused product.

## SECTION 14 - TRANSPORT INFORMATION

<b>UN Number</b>	2222
<b>UN Proper shipping name</b>	ANISOLE
<b>Transport hazard class(es)</b>	3 Flammable liquids
<b>Packing group</b>	III
<b>Limited quantity index</b>	5L
<b>ERAP Index</b>	-
<b>Special precautions</b>	-

## SECTION 15 - REGULATORY INFORMATION

<b>WHIMS CANADA</b>	Skin corrosion/irritation - Skin irritation category 2 Flammable liquids category 3 Specific target organ toxicity - Single exposure category 3 Serious eye damage/eye irritation - Eye irritation category 2
---------------------	--

## SECTION 16 - OTHER INFORMATION

### Further information

The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. It does not represent any guarantee of the properties of the product. Laboratoire MAT Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.

Last Update: 3/18/2024