




Centre Anti-Poison pour le Québec: (800) 463-5060
Tél. (Qc): (418) 660-8666 / 800-890-8666
Fax. (Qc): (418) 660-8998

SAFETY DATA SHEET

SECTION 01 - PRODUCT AND COMPANY IDENTIFICATION

Product Identifier STÉARATE DE CALCIUM		Product Use Laboratory use	
Chemical formula $\text{Ca}[\text{CH}_3(\text{CH}_2)_{16}\text{CO}_2]_2$		Product code CP-0823	Molar weight 607,04
Chemical name / Commercial name / Synonymous Calcium distearate; Sel calcique de l'acide octadécanoïque; Sel calcique de l'acide stéarique; Calcium octadecanoate; Octadecanoic Acid, Calcium Salt; Stearic acid, calcium salt			
Supplier's name Laboratoire MAT		Address-Street 610, Adanac Street	
City Québec		Province Québec	
Postal code G1C 7B7	Internet www.labmat.com	Phone number 418-660-8666 / 800-890-8666	
Emergency phone	CANUTEC: 613-996-6666	CENTRE ANTI-POISON DU QUÉBEC 800-463-5060	
Date SDS 1/15/2020	SDS Prepared by Laboratoire MAT	E-Mail labmat@labmat.com	

SECTION 02 - HAZARDS IDENTIFICATION

Classification WHIMS / GHS	Combustible dusts category 1
Signal Word	WARNING
Hazards statements (H)	May form combustible dust concentrations in air.
Precautionary statements (P)	Keep away from heat, hot surfaces, sparks, open flames and other sources of ignition. No smoking.
Other dangers	NFPA (Risk: 0=No risk; 1=Slight; 2=Moderate; 3=Signifiant; 4=Extreme)
	Health 2 Fire 1 Reactivity 0 Special danger

SECTION 03 - COMPOSITION/INFORMATION ON INGREDIENTS

Ingrédients (Dénomination chimique / synonymes)	Numéro CAS et tout identificateur unique	Concentration (%)
STÉARATE DE CALCIUM	1592-23-0	<=100

SECTION 04 - FIRST AID MEASURES

Eye contact	Wash eyes with large amounts of water for at least 15 minutes while holding eyelids apart to rinse eyes. If irritation persists, seek medical attention.
Skin contact	Wash skin with plenty of water for at least 15 minutes. Remove soiled clothing. If irritation persists, seek medical attention.
Inhalation	If breathing is difficult, give oxygen. If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.
Ingestion	If the person is conscious, give water to drink. Never give anything by mouth to an unconscious person. Consult a physician.
Most important symptoms and effects (acute and delayed)	Ref. section 11.
Immediate medical attention and special treatment, if necessary	In case of medical consultation, keep this sheet available.
General advice	Show this safety data sheet to the doctor in attendance.

SECTION 05 - FIREFIGHTING MEASURES

Flammability	No
Ignition conditions	Fine dust in sufficient concentration may be combustible, or explode if confined to a small space and subject to a source of ignition.
Suitable extinguishing media	Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.
Unsuitable extinguishing media	Data not available.
Hazardous combustion / decomposition products	Hazardous decomposition products formed under fire conditions. - Calcium oxide. Carbon oxides.
Special fire and explosion hazards	May form concentrations combustible dust in the air. May react violently with incompatible products (Ref Section 10).
Special protective equipment and precautions for firefighters	Discard incompatible substances if this can be done without risk. Firefighters should be equipped with standard protective equipment, fireproof clothing, face mask, gloves, protective boots and, where appropriate, self-contained breathing apparatus.

SECTION 06 - ACCIDENTAL RELEASE MEASURES

Methods and materials for containment and cleaning up / Personnel precautions, protective equipment	Evacuate personnel to safe areas. Cut off all sources of ignition. Ensure adequate ventilation. Use personal protective equipment. Pick up with a shovel or broom, taking care not to scatter dust. Avoid dust formation. Clean and rinse with water. Dispose of residues in a container provided for the disposal of hazardous materials. Do not let product enter drains.
--	---

SECTION 07 - HANDLING AND STORAGE

Conditions for safe storage	Keep container tightly closed in a dry and well-ventilated place.
Methods of handling	Avoid breathing dust, fumes, gases, mists, vapors, aerosols. Avoid dust formation. Ensure good ventilation. Wear personal protective equipment. Avoid ingestion and inhalation.

SECTION 08 - EXPOSURE CONTROLS/PERSONAL PROTECTION

Workplace control parameters

Components	CAS-No.	Value	Control parameters	Basis
Calcium stearate	1592-23-0	TWA	10.000000 mg/m3	Canada. Alberta, Occupational Health and Safety Code (table 2: OEL)
Remarks	Occupational exposure limit is based on irritation effects and its adjustment to compensate for unusual work schedules is not required			
		TWAEV	10.000000 mg/m3	Canada. Ontario OELs
		TWA	10.000000 mg/m3	Canada. British Columbia OEL
Does not include stearates of toxic metals.				

Data source	Sigma-Aldrich.
Ventilation	Fan.
Respiratory	If the permissible levels are exceeded, use a mechanical filter / cartridge against NIOSH vapors or a respirator with air supply.
Gloves	Handle with gloves.
Eyes	Safety goggles with safety shutters.
Shoes	Safety shoes.
Clothing	Labcoat. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.
Engineering control	Have safety showers and eyewash stations in the workplace in case of an emergency and a ventilation system to maintain the level of concentrations in the air below the exposure limit values.

SECTION 09 - PHYSICAL AND CHEMICAL PROPERTIES

Physical state	Solid.
Appearance	Poudre ou granules, blanc à jaunâtre.
Odour	Faible odeur d'acide gras.
Odour threshold	Data not available
pH	Donnée non disponible.
Melting point / Freezing point	145-160°C
Initial boiling point	Data not available
Boiling range	Data not available
Flash point	Data not available
Evaporation rate	Data not available
Flammability	No
Lower flammable / Explosive limit	Data not available
Upper flammable / Explosive limit	Data not available
Vapour pressure	Data not available
Vapour density	Data not available
Relative density	Data not available
Solubility	Insoluble dans l'eau: 0.04g/L..
Partition coefficient water/n-octanol	Data not available
Auto-ignition temperature	399°C
Decomposition temperature	Data not available
Viscosity	Data not available

SECTION 10 - STABILITY AND REACTIVITY

Reactivity	Non-reactive under normal conditions.
Chemical stability	Stable under recommended storage conditions.
Possibility of hazardous reactions	May form combustible dust concentrations in the air.
Conditions of instability (Including sensitivity to shock / static discharge / vibration)	Avoid dust formation. Avoid contact with incompatible materials.
Incompatible material	Strong oxidizing agents, strong acids. Strong bases. Alloys based on copper, iron or nickel.
Hazardous decomposition products	Toxic vapors of hydrogen chloride. Carbon oxides. - Calcium oxide.

SECTION 11 - TOXICOLOGICAL INFORMATION

Routes of exposure	Dermal contact.
Acute exposition effects / symptoms:	By exposure route below.
- Eyes	May cause eye irritation. To our knowledge, the product has not been fully studied.
- Skin	May cause skin irritation. To our knowledge, the product has not been fully studied.
- Inhalation	May cause respiratory tract irritation. To our knowledge, the product has not been fully studied.
Acute toxicity (Ingestion)	No known health effects. To our knowledge, the product has not been fully studied.
Chronic exposure effects / symptoms	Animal: skin contact (emulsion): decrease in body weight; intratracheal administration (emulsion): possibility of lung damage and damage to blood vessels in lung tissue.
DL50 (specify species and route of entry)	LD50 Oral - Rat - > 2 000 mg/kg. LD50 Dermal - Rat - 2000mg/kg.
CL50 (specify species and route of entry)	LC50 - Inhalation - Data not available.

SECTION 12 - ECOLOGICAL INFORMATION

Ecotoxicity	Not considered to be harmful to aquatic life.
Persistence and degradability	Readily biodegradable (based on similar product(s)).
Bioaccumulative potential	Does not bioaccumulate.
Mobility in soil	Data not available.
Other adverse effects	Data not available.

SECTION 13 - DISPOSAL CONSIDERATIONS

Waste Disposal Method	Dispose of contents / container in accordance with local / regional / national / international regulations / or contact a specialist waste disposal company.
Contaminated Packaging	Dispose of as unused product.

SECTION 14 - TRANSPORT INFORMATION

UN Number	N/R
UN Proper shipping name	
Transport hazard class(es)	
Packing group	
Limited quantity index	
ERAP Index	
Special precautions	

SECTION 15 - REGULATORY INFORMATION

WHIMS CANADA	Combustible dusts category 1
--------------	------------------------------

SECTION 16 - OTHER INFORMATION

Further information

The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. It does not represent any guarantee of the properties of the product. Laboratoire MAT Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.

Last Update: 1/15/2020