



SAFETY DATA SHEET

SECTION 01 - IDENTIFICATION

Product Identifier	EPINEPHRINE (ADRENALINE)
Other identification (Chemical name, Commercial name, Synonymous)	epinephrine; adrenaline; l-Adrenaline; Adrenalin; Levoepinephrine.
Product code	EP-0113
Chemical formula	$C_9H_{13}NO_3$
Molar weight	183.21
Recommended use and Restrictions on use	For laboratory, school, commercial or industrial use. Not for medical or household use. Do not use for medical, food or household purposes.
Supplier	LABORATOIRE MAT 610, rue Adanac Québec Québec G1C 7B7 418-660-8666 Mon-Fri 8h-16h www.labmat.com labmat@labmat.com
Emergency phone	418-660-8666 Mon-Fri 8h-16h CENTRE ANTI-POISON DU QUÉBEC 800-463-5060
Date SDS	2025-06-17

SECTION 02 - HAZARDS IDENTIFICATION

WHIMS CANADA

- Acute toxicity - Dermal - category 2
- Acute toxicity - Inhalation - category 3
- Acute toxicity - Oral - category 3

PICTOGRAMS



Signal Word

DANGER

Hazards statements (H)

- Fatal in contact with skin
- Toxic if inhaled
- Toxic if swallowed

Precautionary statements (P)

- Do not get in eyes, on skin, or on clothing.
- Wear protective gloves (nitrile, butyl, neoprene), protective clothing and eye and face protection.
- IF ON SKIN: Wash with plenty of water
- Immediately call a POISON CENTER or a physician.
- Store locked up.
- Dispose of contents and container in accordance with local, regional and national regulations, or contact a specialist waste disposal company.
- Avoid breathing mists, gases, vapors and other fumes, or the product itself.
- Use only outdoors or in a well-ventilated area.
- IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing.
- Call a POISON CENTER or a physician.
- Store in a well ventilated place. Keep container tightly closed.
- Rinse mouth.

Other dangers

NFPA (Risk: 0=No risk; 1=Slight; 2=Moderate; 3=Significant; 4=Extreme)

Health 4
Fire 0
Reactivity 0
Special danger

SECTION 03 - COMPOSITION/INFORMATION ON INGREDIENTS

Component	No. CAS	% Weight
EPINEPHRINE (ADRENALINE)	51-43-4	<=100%

SECTION 04 - FIRST AID MEASURE

Eye contact	Wash eyes with large amounts of water for at least 15 minutes while holding eyelids apart to rinse eyes. If irritation persists, seek medical attention.
Skin contact	Wash with soap using a large amount of water. Immediately transport the victim to the hospital. If irritation persists, seek medical attention.
Inhalation	Move the unwell person to the fresh air. If breathing is difficult, give oxygen. Consult a physician.
Ingestion	If the person is conscious, give water to drink. Never give anything by mouth to an unconscious person. Consult a physician.
Most important symptoms and effects (acute and delayed)	Main symptoms of high exposure: Headaches. Dizziness. Dizziness. Fatigue. Breathing difficulties. Tremors. Neurological damage. Ref. section 11.
Immediate medical attention and special treatment, if necessary	Treat according to symptoms. Show this sheet to the attending physician.

SECTION 05 - FIREFIGHTING MEASURES

Suitable extinguishing media	Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.
Unsuitable extinguishing media	Data not available.
Combustion products	Hazardous combustion products formed under fire conditions: Carbon oxides. Nitrogen oxides (NOx).
Specific hazards of the dangerous product	May react violently with incompatible products (Ref Section 10).
Special protective equipment and precautions for firefighters	Discard incompatible substances if this can be done without risk. Firefighters should be equipped with standard protective equipment, fireproof clothing, face mask, gloves, protective boots and, where appropriate, self-contained breathing apparatus.

SECTION 06 - ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency measures	Evacuate personnel to safe areas. When handling, wear appropriate safety equipment (reference Section 8 for protective equipment to be used). Ensure a good ventilation. Use NIOSH cartridge respiratory protection if necessary or for larger spills. Avoid dust formation. Avoid breathing in vapours, spray mists or gases.
Methods and materials for containment and cleaning up	Pick up with a shovel or broom, taking care not to scatter dust. Dispose of residues in a container for disposal of hazardous materials. When handling, wear suitable safety equipment. Discharge into the environment must be avoided.

SECTION 07 - HANDLING AND STORAGE

Conditions for safe storage	Recommended storage temperature: 2-8 ° C. Keep container tightly closed in a dry, well-ventilated place. Protect from light and sunlight. Store under inert gas. Keep container tightly closed and store away from heat, moisture, and incompatible products (ref. section 10).
Methods of handling	Avoid contact with eyes, skin and clothing. Avoid ingestion and inhalation. Wear personal protective equipment (ref. section 8) when handling. Always ensure good ventilation. Apply the usual standard hygiene rules: Wash your hands after use. Do not eat or drink during use. Avoid the formation of dust.

SECTION 08 - EXPOSURE CONTROLS/PERSONAL PROTECTION

WORKPLACE CONTROL PARAMETERS

Component	CAS-No.	Value
EPINEPHRINE (ADRENALINE)	51-43-4	No occupational exposure limits established by region-specific regulators - Quebec, Alberta, Ontario, British Columbia.

Respiratory	If work under the hood is not possible, or if the permissible levels are exceeded, use NIOSH cartridge respiratory protection, or an air-supplied respirator.
Gloves	Handle with protective gloves. Suggested material: Nitrile. Butyl. Neoprene. The type, thickness and length of the glove must be chosen according to the use, the concentration of the product, as well as the duration of use. Replace gloves regularly for better protection.
Eyes	Safety goggles with safety shutters.
Shoes	Use safety shoes.
Clothes	Labcoat. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.
Engineering control	Use fan. Have safety showers and eyewash stations in the workplace in case of an emergency and a ventilation system to maintain the level of concentrations in the air below the exposure limit values.

SECTION 09 - PHYSICAL AND CHEMICAL PROPERTIES

Physical state	Solid.
Color	Beige.
Odour	Data not available.
Odour threshold	Data not available.
Melting point and freezing point	215 °C (419 °F).
Boiling point	Data not available.
Flammability	Data not available.
Lower flammable / Explosive limit	Data not available.
Upper flammable / Explosive limit	Data not available.
Flash point	Data not available.
Auto-ignition temperature	Data not available.
Decomposition temperature	Data not available.
pH	Data not available.
Kinematic viscosity	Data not available.
Solubility	Data not available.
Partition coefficient water/n-octanol	Data not available.
Vapour pressure	Data not available.
Relative density	Data not available.
Vapour density	Data not available.
Particle characteristics	Crystalline powder.

SECTION 10 - STABILITY AND REACTIVITY

Reactivity	May react violently with incompatible substances.
Chemical stability	Stable under recommended storage conditions.
Possibility of hazardous reactions	May react violently with incompatible substances.
Conditions to avoid	Avoid moisture. Avoid contact with incompatible materials and extreme temperatures.
Incompatible materials	Strong acids. Acid anhydrides. Acid chlorides. Oxydants. Bases. Copper. Iron. Zinc. Metals.
Hazardous decomposition products	Carbon oxides. Nitrogen oxides (NO _x).

SECTION 11 - TOXICOLOGICAL INFORMATION

EPINEPHRINE (ADRENALINE)

Routes of exposure	Ingestion, inhalation, skin and eye contact.
Acute exposition effects / symptoms:	By exposure route below.
- Eyes	Irritation.
- Skin	Can be fatal in case of absorption through the skin. Irritation.
- Inhalation	Irritation of the mucous membranes and respiratory tract.
Acute toxicity (Ingestion)	To our knowledge, the chemical, physical and toxicological properties have not been fully investigated.
Chronic exposure effects / symptoms	Headache. Dizziness. Dizziness. Fatigue. Muscle weakness. Loss of consciousness. Nervous system damage. Conjunctivitis. Sweating. Liver and kidney damage. Confusion. Convulsions. Fever.
DL50 (specify species and route of entry)	LD50 Oral - Data not available. LD50 Dermal - Rat - 62 mg/kg.
CL50 (specify species and route of entry)	LC50 Inhalation - Data not available.

SECTION 12 - ECOLOGICAL INFORMATION

EPINEPHRINE (ADRENALINE)

Ecotoxicity	Data not available.
Persistence and degradability	Data not available.
Bioaccumulative potential	Data not available.
Mobility in soil	Data not available.
Other adverse effects	Data not available.

SECTION 13 - DISPOSAL CONSIDERATIONS

Waste Disposal Method	Dispose of contents and container in accordance with local, regional and national regulations, or contact a specialist waste disposal company.
Contaminated Packaging	Dispose of as unused product.

SECTION 14 - TRANSPORT INFORMATION

UN Number	2811
UN Proper shipping name	SOLIDE ORGANIQUE TOXIQUE, N.S.A. (Epinephrine (adrenaline))
Transport hazard class(es)	Matières toxiques 6.1
Packing group	II
Limited quantity index	0,5 kg
ERAP Index	-
Special precautions	16

SECTION 15 - REGULATORY INFORMATION

WHIMS CANADA	- Acute toxicity - Dermal - category 2 - Acute toxicity - Inhalation - category 3 - Acute toxicity - Oral - category 3
--------------	--

SECTION 16 - OTHER INFORMATION

CNESST: Commission des normes, de l'équité et de la santé et sécurité au travail
NIH: National institute of health (U.S. National Library of Medicine)
ECHA: Agence Européenne de Chimie
ACGIH : American Conference of Governmental Industrial Hygienists
AIHA : American Industrial Hygiene Association
VECD: Valeur d'exposition courte durée
VEMP: Valeur d'exposition moyenne pondérée
CMRG : Chemical Manufacturer's Recommended Guidelines
OSHA : United States Department of Labor - Occupational Safety and Health Administration
TLV : Threshold limit value
TWA: Time-Weighted-Average
STEL: Short Term Exposure Limit
CEIL: Ceiling
RSST: Règlement sur la santé et sécurité au travail (Québec)
INRS: l'Institut national de recherche et de sécurité pour la prévention des accidents du travail et des maladies professionnelles (France)

The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. It does not represent any guarantee of the properties of the product. Laboratoire MAT Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.

LAST UPDATE: 2025-06-17