



Centre Anti-Poison pour le Québec: (800) 463-5060

Tél. (Qc): (418) 660-8666 / 800-890-8666


Fax. (Qc): (418) 660-8998

SAFETY DATA SHEET

SECTION 01 - PRODUCT AND COMPANY IDENTIFICATION

Product Identifier ERIOCHROME® NOIR T 0.5%w/v in TRIÉTHANOLAMINE		Product Use Laboratory use	
Chemical formula -		Product code ES-0105	Molar weight
Chemical name / Commercial name / Synonymous TRIETHANOLAMINE, 2,2',2''-NITRILOTRIETHANOL, TRIHYDROXYTRIETHYLAMINE, TRIS (HYDROXYETHYL) AMINE, TRIETHYLOLAMINE, TROLAMINE, DALTOGEN, STEROLAMIDE, THIOFACO T-35			
Supplier's name Laboratoire MAT		Address-Street 610, Adanac Street	
City Québec		Province Québec	
Postal code G1C 7B7	Internet www.labmat.com	Phone number 418-660-8666 / 800-890-8666	
Emergency phone	CANUTEC: 613-996-6666 CENTRE ANTI-POISON DU QUÉBEC 800-463-5060		
Date SDS 10/12/2022	SDS Prepared by Laboratoire MAT	E-Mail labmat@labmat.com	

SECTION 02 - HAZARDS IDENTIFICATION

Classification WHIMS / GHS	Skin corrosion/irritation - Skin irritation category 2 Respiratory or skin sensitization - Skin sensitize category 1B Specific target organ toxicity - Single exposure category 3 Serious eye damage/eye irritation - Eye irritation category 2	
Signal Word	WARNING	
Hazards statements (H)	H315 Causes skin irritation. H317 May cause an allergic skin reaction. H319 Causes serious eye irritation. H335 May cause respiratory irritation.	
Precautionary statements (P)	P261 Avoid breathing dust / fume / gas / mist / vapors / spray. P264 Wash the areas of the body that have been in contact with the product after handling. P271 Use only outdoors or in a well-ventilated area. P272 Contaminated work clothing should not be allowed out of the workplace. P280 Wear protective gloves/protective clothing/eye protection/face protection. P302 + P352 IF ON SKIN: Wash with plenty of soap and water. P304 + P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P312 Call a POISON CENTER or doctor/physician if you feel unwell. P321 Specific treatment (see section 4 of the SDS and on this label). P332 + P313 If skin irritation occurs: Get medical advice/attention. P333 + P313 If skin irritation or rash occurs: Get medical advice/attention. P337 + P313 If eye irritation persists: Get medical advice/attention. P362 + P364 Take off contaminated clothing and wash it before reuse. P403 + P233 Store in a well-ventilated place. Keep container tightly closed. P405 Store locked up. P501 Dispose of contents/container in accordance with local / regional / national / international regulations or contact a specialist waste disposal company.	
PICTOGRAMS		
Other dangers	NFPA (Risk: 0=No risk; 1=Slight; 2=Moderate; 3=Signifiant; 4=Extreme)	
	Health 2 Fire 1 Reactivity 1 Special danger	

SECTION 03 - COMPOSITION/INFORMATION ON INGREDIENTS

Ingrédients (Dénomination chimique / synonymes)	Numéro CAS et tout identificateur unique	Concentration (%)
Triéthanolamine	102-71-6	99.6
Eriochrome Noir T	1787-61-7	0.4

SECTION 04 - FIRST AID MEASURES

Eye contact	Wash eyes with large amounts of water for at least 15 minutes while holding eyelids apart to rinse eyes. If irritation persists, seek medical attention.
Skin contact	Wash skin with plenty of water for at least 15 minutes. Remove soiled clothing. If irritation persists, seek medical attention.
Inhalation	Move the unwell person to the fresh air. If breathing is difficult, give oxygen. Consult a physician.
Ingestion	Do NOT induce vomiting. If the person is conscious, rinse the mouth with water. Never give anything by mouth to an unconscious person. Consult a physician.
Most important symptoms and effects (acute and delayed)	Main symptoms of high exposure: Headaches. Skin, eye and respiratory system irritation. Vomiting. Redness. Ref. section 11.
Immediate medical attention and special treatment, if necessary	In case of medical consultation, keep this sheet available.
General advice	Show this safety data sheet to the doctor in attendance.

SECTION 05 - FIREFIGHTING MEASURES

Flammability	No
Ignition conditions	Non flammable. But remains combustible. Strong oxidizing agents, heat, sparks and open flame.
Suitable extinguishing media	Use water spray, alcohol-resistant foam, or dry chemical.
Unsuitable extinguishing media	Data not available.
Hazardous combustion products	Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NO _x), Sulphur oxides, Sodium oxides. Hydrogen cyanide (hydrocyanic acid). Formaldehyde.
Special fire and explosion hazards	May react violently with incompatible products (Ref Section 10).
Special protective equipment and precautions for firefighters	Discard incompatible substances if this can be done without risk. Firefighters should be equipped with standard protective equipment, fireproof clothing, face mask, gloves, protective boots and, where appropriate, self-contained breathing apparatus.

SECTION 06 - ACCIDENTAL RELEASE MEASURES

Methods and materials for containment and cleaning up / Personnel precautions, protective equipment	Evacuate personnel to safe areas. Absorb the product with sand or vermiculite. Dilute residues with water, clean and rinse. Ensure a good ventilation of the premises. Dispose of residues in a container for disposal of hazardous materials. When handling, wear suitable safety equipment. Use breathing apparatus if necessary. Avoid breathing vapors, mist or gas. Discharge into the environment must be avoided.
--	--

SECTION 07 - HANDLING AND STORAGE

Conditions for safe storage	Hygroscopic. Store in a well-ventilated area. Store in a cool, dry place. Keep container tightly closed and store away from heat, moisture, and incompatible products. Protect from light and sunlight.
Methods of handling	Keep away from moisture. Avoid ingestion and inhalation. Provide appropriate exhaust ventilation at places where dust or vapor is formed. Avoid contact with the skin, eyes and clothes. Wear personal protective equipment when handling. Always ensure good ventilation. Transport according to TDG (ref Section 14)

SECTION 08 - EXPOSURE CONTROLS/PERSONAL PROTECTION

Workplace control parameters

Composants	No.-CAS	Valeur	Paramètres de contrôle	Base
2,2',2''-Nitrilotriethanol	102-71-6	TWAEV	0.500000 ppm 3.100000 mg/m3	Canada. Ontario OELs
		TWA	5.000000 mg/m3	Canada. LEP Colombie Britannique
		TWA	5.000000 mg/m3	Canada. Alberta, Code de santé et de sécurité au travail (tableau 2 : VLE)
Remarques	La limite d'exposition professionnelle est basée sur les effets de l'irritation et son ajustement pour compenser les emplois du temps de travail inhabituels n'est pas nécessaire			
		VEMP	5.000000 mg/m3	Québec. Règlement sur la santé et la sécurité du travail, Annexe 1 Partie 1: Valeurs d'exposition admissibles des contaminants de l'air
	Sensibilisant			
	VEMP			
	5 mg/m3 Québec. Règlement sur la santé et la sécurité du travail, Annexe 1 Partie 1: Valeurs d'exposition admissibles des contaminants de l'air			
	Sensibilisant			
		LMPT	0.500000 ppm 3.100000 mg/m3	Tableau de l'Ontario: Limites d'exposition professionnelle pris en vertu de la loi sur la santé et la sécurité au travail.
		LMPT	0.5 ppm 3.1 mg/m3	Tableau de l'Ontario: Limites d'exposition professionnelle pris en vertu de la loi sur la santé et la sécurité au travail.
Components	CAS-No.	Control	Value	Basis
Eriochrome Black T (Indicator)	1787-61-7	TLV, TWA, STEL	No data available	Canada. Alberta, Occupational Health and Safety Code (table 2: OEL)
		TLV, TWA, STEL	No data available	Canada. British Columbia OEL
		TLV, TWA, STEL	No data available	Québec. Regulation respecting occupational health and safety, Schedule 1, Part 1: Permissible exposure values for airborne contaminants

Data source	Sigma-Aldrich.
Ventilation	Fan.
Respiratory	If the permissible levels are exceeded, use a mechanical filter / cartridge against NIOSH vapors or a respirator with air supply.
Gloves	Handle with gloves.
Eyes	Safety goggles with safety shutters.
Shoes	Safety shoes.
Clothing	Labcoat. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.
Engineering control	Have safety showers and eyewash stations in the workplace in case of an emergency and a ventilation system to maintain the level of concentrations in the air below the exposure limit values.

SECTION 09 - PHYSICAL AND CHEMICAL PROPERTIES

Physical state	Liquid.
Appearance	Liquide visqueux bleu foncé - noir.
Odour	Odeur légèrement d'ammoniac.
Odour threshold	Data not available
pH	Solution aqueuse Triéthanolamine 0.1N = pH 10.5.
Melting point / Freezing point	17.9-21°C (Triéthanolamine).
Initial boiling point	335°C (Triéthanolamine).
Boiling range	Data not available
Flash point	179 - 193 °C @ 101.3 - 101.325 kPa (Triéthanolamine).
Evaporation rate	Data not available
Flammability	No
Lower flammable / Explosive limit	1.3% (Triéthanolamine).
Upper flammable / Explosive limit	8.5% (Triéthanolamine).
Vapour pressure	0.01 mmHg @ 20°C (Triéthanolamine).
Solubility	Miscible avec l'eau en toutes proportions.
Vapour density	5.15 (Air 1.0) (Triéthanolamine)-
Relative density	1.124 g/ml à 25 °C.
Partition coefficient water/n-octanol	log Kow = -1.00 (Triéthanolamine)-
Auto-ignition temperature	325°C (Triéthanolamine).
Decomposition temperature	Data not available
Viscosity	600 mPa.s à 25°C (Triéthanolamine).

SECTION 10 - STABILITY AND REACTIVITY

Reactivity	Basic product, reacts violently with strong acids.
Chemical stability	Stable under recommended storage conditions. Hygroscopic.
Possibility of hazardous reactions	May react violently with incompatible substances.
Conditions of instability (Including sensitivity to shock / static discharge / vibration)	Avoid exposure to air, light and moisture. Avoid contact with incompatible materials. Avoid excessive heat.
Incompatible material	Halogenated organic compounds. Strong oxidizing agents, metals. Isocyanates. Acids. Peroxides. Halogens. Anhydrides of acid. Strong reducing agents. Hydrides. The acid halides.
Hazardous decomposition products	Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NOx), Sulphur oxides, Sodium oxides. Hydrogen cyanide (hydrocyanic acid). Formaldehyde.

SECTION 11 - TOXICOLOGICAL INFORMATION

TRIETHANOLAMINE

Routes of exposure	Ingestion, inhalation, skin and eyes.
Acute exposition effects / symptoms:	By exposure route below.
- Eyes	Severe irritation.
- Skin	Irritation and dermatitis.
- Inhalation	Irritation of the mucous membranes and respiratory tract.
Acute toxicity (Ingestion)	Irritation of the mucous membranes. Dysphagia, abdominal pain, cramps, diarrhea, headache, dizziness, diaphoresis, convulsions, nausea and vomiting.
Chronic exposure effects / symptoms	Burning sensation, dermatitis, conjunctivitis, nervous disorders, liver and kidney damage, chest pain, cough, dyspnea, laryngitis, headache, dizziness, diaphoresis, asthenia, weight loss and loss of appetite, nausea and vomiting. Is recognized as an unclassifiable agent for its carcinogenicity (Group 3 according to IARC).
DL50 (specify species and route of entry)	LD50 Oral - Rat - 6400 mg/kg. LD50 Dermal - Rat - > 16 mL/kg LD50 Dermal - Rabbit - 2000 mg/kg.
CL50 (specify species and route of entry)	LC50 - Inhalation - Data not available.

ERIOCHROME® BLACK T (INDICATOR)

Routes of exposure	Ingestion, inhalation, skin and eyes.
Acute exposition effects / symptoms:	By exposure route below.
- Eyes	Irritation.
- Skin	Irritation.
- Inhalation	Irritation of the mucous membranes and respiratory tract. Nervous disorders, dyspnea, cough, headache, nausea and vomiting.
Acute toxicity (Ingestion)	Irritation of the mucous membranes. Gastrointestinal disorders, cramps, diarrhea, headache, nausea and vomiting.
Chronic exposure effects / symptoms	Burning sensation, nervous disorders, cough, dyspnea, laryngitis, headache, fatigue, loss of appetite, nausea and vomiting.
DL50 (specify species and route of entry)	LD50 Oral - Rat - 17,590 mg/kg Remarks: Behavioral:Excitement. Behavioral:Change in motor activity (specific assay). LD50 Dermal - Rabbit - > 3000 mg / kg
CL50 (specify species and route of entry)	LC50 Inhalation - Rat - 4h - > 2 mg / L

SUMMARY

Acute exposure effects / Symptoms:	By exposure routes below.
Ingestion	To our knowledge, the product has not been fully evaluated
Inhalation	To our knowledge, the product has not been fully evaluated
Skin	To our knowledge, the product has not been fully evaluated
Eyes	To our knowledge, the product has not been fully evaluated
Chronic exposure effects / Symptoms:	To our knowledge, the product has not been fully evaluated
ETA Mix (Estimated Acute Toxicity)	LD50 Oral: 6418mg/kg - Rat LD50 Dermal: 2003 mg/kg - Rabbit LC50 Inhalation: No data available

SECTION 12 - ECOLOGICAL INFORMATION

Ecotoxicity	Triethanolamine: EC50 - Daphnia (water flea) - 609.98 mg/l - 48 h LC50 - Pimephales promelas (fathead minnow) - 10600 - 13000 mg/L - 96 h
Persistence and degradability	Soluble in water. Persistence is unlikely based on information available.
Bioaccumulative potential	Data not available.
Mobility in soil	Probable mobility in the environment due to its solubility in water.
Other adverse effects	Data not available.

SECTION 13 - DISPOSAL CONSIDERATIONS

Waste Disposal Method	Dispose of contents / container in accordance with local / regional / national / international regulations / or contact a specialist waste disposal company.
Contaminated Packaging	Dispose of as unused product.

SECTION 14 - TRANSPORT INFORMATION

UN Number	N/R
UN Proper shipping name	
Transport hazard class(es)	
Packing group	
Limited quantity index	
ERAP Index	
Special precautions	

SECTION 15 - REGULATORY INFORMATION

WHIMS CANADA	Skin corrosion/irritation - Skin irritation category 2 Respiratory or skin sensitization - Skin sensitize category 1B Specific target organ toxicity - Single exposure category 3 Serious eye damage/eye irritation - Eye irritation category 2
--------------	--

SECTION 16 - OTHER INFORMATION

Further information

The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. It does not represent any guarantee of the properties of the product. Laboratoire MAT Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.

Last Update: 10/12/2022