




SAFETY DATA SHEET

SECTION 01 - PRODUCT AND COMPANY IDENTIFICATION

Product Identifier FORMALDEHYDE 0.02N		Product Use Laboratory use	
Chemical formula HCHO		Product code FS-2227	Molar weight 30,03
Chemical name / Commercial name / Synonymous FORMALDEHYDE, ALDÉHYDE FORMIQUE, ALDÉHYDE DE MÉTHYLE, OXYDE DE MÉTHYLÈNE, OXOMÉTHANE, OXYMÉTHYLÈNE, POLYOXYMÉTHYLÈNE GLYCOL, METHANAL, MORBICID			
Supplier's name Laboratoire MAT		Address-Street 610, Adanac Street	
City Québec		Province Québec	
Postal code G1C 7B7	Internet www.labmat.com	Phone number 418-660-8666 Lun-Ven 8h-16h	
Emergency phone	418-660-8666 Lun-Ven 8h-16h	CENTRE ANTI-POISON DU QUÉBEC 800-463-5060	
Date SDS 1/27/2025	SDS Prepared by Laboratoire MAT	E-Mail labmat@labmat.com	

SECTION 02 - HAZARDS IDENTIFICATION

Classification WHIMS / GHS	Not a hazardous substance according to WHMIS 2015
Other dangers	NFPA (Risk: 0=No risk; 1=Slight; 2=Moderate; 3=Signifiant; 4=Extreme)
	Health 2 Fire 0 Reactivity 0 Special danger

SECTION 03 - COMPOSITION/INFORMATION ON INGREDIENTS

Ingrédients (Dénomination chimique / synonymes)	Numéro CAS et tout identificateur unique	Concentration (%P/P)
Formaldéhyde	50-00-0	0.06
Méthanol (comme stabilisant)	67-56-1	0.01

SECTION 04 - FIRST AID MEASURES

Eye contact	Wash eyes with large amounts of water for at least 15 minutes while holding eyelids apart to rinse eyes. If irritation persists, seek medical attention.
Skin contact	Wash skin with plenty of water for at least 15 minutes. Remove soiled clothing. If irritation persists, seek medical attention.
Inhalation	If breathed in, move person into fresh air. If breathing is difficult, give oxygen. Consult a physician.
Ingestion	If the person is conscious, give water to drink. Never give anything by mouth to an unconscious person. Consult a physician.
Most important symptoms and effects (acute and delayed)	To our knowledge, the chemical, physical and toxicological properties have not been fully investigated. Ref. section 11.
Immediate medical attention and special treatment, if necessary	Treat according to symptoms. Show this sheet to the attending physician.

SECTION 05 - FIREFIGHTING MEASURES

Flammability	No
Suitable extinguishing media	Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.
Unsuitable extinguishing media	Data not available.
Hazardous combustion products	Hazardous combustion products formed under fire conditions: Carbon oxides.
Specific hazards of the dangerous product	To our knowledge, the product has not been fully evaluated. When concentrated, the product reacts according to the following characteristics: Combustible liquid. Moderate risks of fire or explosion in the presence of heat or flame. Vapors may form flammable or explosive mixtures with air. Violent reactions with perchloric acid + aniline, performic acid, magnesium carbonate, nitromethane and hydrogen peroxide. May react violently with incompatible products (Ref Section 10).
Special protective equipment and precautions for firefighters	Discard incompatible substances if this can be done without risk. Firefighters should be equipped with standard protective equipment, fireproof clothing, face mask, gloves, protective boots and, where appropriate, self-contained breathing apparatus.

SECTION 06 - ACCIDENTAL RELEASE MEASURES

Methods and materials for containment and cleaning up / Personal precautions, protective equipment	Evacuate personnel to safe areas. Ensure a good ventilation of the premises. When handling, wear appropriate safety equipment. Use NIOSH cartridge respiratory protection for larger spills. (Reference Section 8 for protective equipment to be used.) Absorb the product with sand or vermiculite. Dilute residues with water, clean and rinse. Dispose of residues in a container provided for the disposal of hazardous materials. Do not let product enter drains.
---	---

SECTION 07 - HANDLING AND STORAGE

Conditions for safe storage	Store in a cool, dry place. Protect from the sun's rays. Keep container tightly closed and store away from heat, water, moisture, and incompatible products (ref. section 10). Use venting and electrical equipment that is grounded and does not produce ignition sources (sparks).
Methods of handling	Always open containers slowly to allow any excess pressure to vent. Avoid ingestion and inhalation. Keep away from sources of ignition - No smoking. Wear personal protective equipment when handling. Always ensure good ventilation. Apply the usual standard hygiene rules: Wash your hands after use. Do not eat or drink during use.

SECTION 08 - EXPOSURE CONTROLS/PERSONAL PROTECTION

Workplace control parameters

Composants	No.-CAS	Control parameters	Value	Basis
Formaldéhyde	50-00-0	(c)	1.000000 ppm 1.300000 mg/m3	Canada. Alberta, Occupational Health and Safety Code (table 2: OEL)
Remarques	cancérogène classe A1 pour l'homme			
		TWA	0.750000 ppm 0.900000 mg/m3	Canada. Alberta, Occupational Health and Safety Code (table 2: OEL)
	cancérogène classe A1 pour l'homme			
		TWA	0.300000 ppm	Canada. British Columbia OEL
	cancérogène classe A1 pour l'homme			
		C	1.000000 ppm	Canada. British Columbia OEL
	cancérogène classe A1 pour l'homme			
		P	1,5 ppm 3.000000 mg/m3	Québec. Règlement sur la santé et la sécurité du travail, Annexe 1 Partie 1: Valeurs d'exposition admissibles des contaminants de l'air
	Cancérogène classe A1 pour l'homme. Une substance dont la recirculation est prohibée conformément à l'article 108.			

Components	CAS-No.	Value	Control parameters	Basis
Methanol	67-56-1	TWA	200.000000 ppm 262.000000 mg/m3	Canada. Alberta, Occupational Health and Safety Code (table 2: OEL)
Remarks	Substance may be readily absorbed through intact skin			
		STEL	250.000000 ppm 328.000000 mg/m3	Canada. Alberta, Occupational Health and Safety Code (table 2: OEL)
	Substance may be readily absorbed through intact skin			
		TWA	200.000000 ppm	Canada. British Columbia OEL
	Contributes significantly to the overall exposure by the skin route.			
		STEL	250.000000 ppm	Canada. British Columbia OEL
	Contributes significantly to the overall exposure by the skin route.			
		TWAEV	200.000000 ppm 262.000000 mg/m3	Québec. Regulation respecting occupational health and safety, Schedule 1, Part 1: Permissible exposure values for airborne contaminants
	Skin (percutaneous)			
		STEV	250.000000 ppm 328.000000 mg/m3	Québec. Regulation respecting occupational

				health and safety, Schedule 1, Part 1: Permissible exposure values for airborne contaminants
	Skin (percutaneous)			
		TWA	200.000000 ppm	USA. ACGIH Threshold Limit Values (TLV)
		STEL	250.000000 ppm	USA. ACGIH Threshold Limit Values (TLV)

Data source	Sigma-Aldrich (Millipore Sigma)
Ventilation	Fan.
Respiratory	If permitted levels are exceeded, use NIOSH cartridge respiratory protection, or an air-supplied respirator.
Gloves	Handle with gloves. Suggested material: Nitrile. Butyle. Neoprene. The type, thickness and length of the glove must be chosen according to the use, the concentration of the product, as well as the duration of use. Replace gloves regularly for better protection.
Eyes	Safety goggles with safety shutters.
Shoes	Safety shoes.
Clothing	Labcoat. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.
Engineering control	Have safety showers and eyewash stations in the workplace in case of an emergency and a ventilation system to maintain the level of concentrations in the air below the exposure limit values.

SECTION 09 - PHYSICAL AND CHEMICAL PROPERTIES

Physical state	Liquid.
Appearance	Liquide clair.
Odour	Forte, piquante..
Odour threshold	Data not available
pH	2.7 - 3.0; 1% @ 25°C.
Melting point / Freezing point	Data not available
Initial boiling point	Data not available
Boiling range	Data not available
Flash point	Data not available
Evaporation rate	Data not available
Flammability	No
Lower flammable / Explosive limit	Data not available
Upper flammable / Explosive limit	Data not available
Vapour pressure	Data not available
Solubility	Miscible avec l'eau, l'alcool et l'acétone.
Vapour density	Data not available
Relative density	1.00g/ml
Partition coefficient water/n-octanol	Data not available
Auto-ignition temperature	Data not available
Decomposition temperature	Data not available
Viscosity	Data not available

SECTION 10 - STABILITY AND REACTIVITY

Reactivity	May react violently with incompatible substances.
Chemical stability	Stable under recommended storage conditions.
Possibility of hazardous reactions	Stable under normal conditions. Vapors may form explosive mixture with air.
Conditions to avoid, including static discharge, shock or vibration	Avoid contact with incompatible materials and extreme temperatures. Heat, flames, sparks. In contact with air, formaldehyde oxidizes to formic acid which is unstabilized and causes polymerization (deposition).
Incompatible materials	When pure, the product reacts with the following products: Acids, Oxidants, Acid Chlorides, Acid Anhydrides, Alkali Metals, Reducing Agents. Strong bases. Isocyanates. Nitriles. Amines. Aniline. Phenol.
Hazardous decomposition products	Hazardous decomposition products formed under fire conditions. Carbon oxides. Paraformaldehyde.

SECTION 11 - TOXICOLOGICAL INFORMATION

FORMALDEHYDE 37%

Routes of exposure	Ingestion, inhalation, skin and eyes.
Acute exposition effects / symptoms:	By exposure route below.
- Eyes	Severe irritation and burns that may cause permanent eye damage.
- Skin	Irritation and dermatitis. Prolonged skin contact may result in an allergic reaction characterized mainly by erythematous or eczematous lesions.
- Inhalation	Irritation of the mucous membranes and respiratory tract. Pains in the chest, respiratory allergies, cough, dyspnea, headache, dizziness, watery eyes, congestion, bronchial spasms and may lead to pulmonary edema.
Acute toxicity (Ingestion)	Irritation and burning of the esophagus and stomach. Abdominal pain, cramps, diarrhea, nausea and vomiting, hematemesis, acidosis, hematuria, anuria, vertigo, pallor, blindness, convulsions, stupor, respiratory collapse, coma and can lead to death.
Chronic exposure effects / symptoms	Is recognized as a carcinogen (class 1) by IARC. Burning sensation, dermatitis, conjunctivitis, chest pain, eye and lung damage, respiratory and skin allergies, cough, dyspnoea, bronchitis, dry throat, headache, dizziness, confusion, irritability, tearing, choking, sleep, intense thirst, sweating, salivation, fatigue, paleness, muscle weakness, weight loss and loss of appetite, convulsions, nausea and vomiting.
DL50 (specify species and route of entry)	LD50 Oral - Rat - 460 mg/kg. LD50 Dermal - Data not available.
CL50 (specify species and route of entry)	LC50 Inhalation - Rat - 4h - 463 ppm.

SUMMARY

Acute exposure effects / Symptoms:	By exposure routes below.
Ingestion	To our knowledge, the product has not been fully evaluated
Inhalation	To our knowledge, the product has not been fully evaluated
Skin	To our knowledge, the product has not been fully evaluated
Eyes	To our knowledge, the product has not been fully evaluated
Chronic exposure effects / Symptoms:	To our knowledge, the product has not been fully evaluated
ETA Mix (Estimated Acute Toxicity)	LD50 Oral: > 5000 mg/kg - Rat LD50 Dermal: - No data available LC50 Inhalation: >100 000 ppm - 4h - Rat

SECTION 12 - ECOLOGICAL INFORMATION

Available ecological information	No
---	----

SECTION 13 - DISPOSAL CONSIDERATIONS

Waste Disposal Method	Dispose of contents and container in accordance with local, regional and national regulations, or contact a specialist waste disposal company.
Contaminated Packaging	Dispose of as unused product.

SECTION 14 - TRANSPORT INFORMATION

UN Number	N/R
UN Proper shipping name	
Transport hazard class(es)	
Packing group	
Limited quantity index	
ERAP Index	
Special precautions	

SECTION 15 - REGULATORY INFORMATION

WHMIS CANADA	Not a hazardous substance according to WHMIS 2015
--------------	---

SECTION 16 - OTHER INFORMATION

Further information

The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. It does not represent any guarantee of the properties of the product. Laboratoire MAT Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.

Last Update: 1/27/2025