



Centre Anti-Poison pour le Québec: (800) 463-5060

Tél. (Qc): (418) 660-8666 / 800-890-8666


Fax. (Qc): (418) 660-8998

SAFETY DATA SHEET

SECTION 01 - PRODUCT AND COMPANY IDENTIFICATION

Product Identifier LITHIUM CHLORIDE		Product Use Laboratory use	
Chemical formula LiCl		Product code LP-0911; LC-0911; LT-0911; LR-0111	Molar weight 42,39
Chemical name / Commercial name / Synonymous CHLORURE DE LITHIUM ANHYDRE, Lithium chloride, Lithium chlorure anhydre, LITHIUM CHLORIDE ANHYDROUS,			
Supplier's name Laboratoire MAT		Address-Street 610, Adanac Street	
City Québec		Province Québec	
Postal code G1C 7B7	Internet www.labmat.com	Phone number 418-660-8666 / 800-890-8666	
Emergency phone	CANUTEC: 613-996-6666 CENTRE ANTI-POISON DU QUÉBEC 800-463-5060		
Date SDS 7/5/2019	SDS Prepared by Laboratoire MAT	E-Mail labmat@labmat.com	

SECTION 02 - HAZARDS IDENTIFICATION

Classification WHIMS / GHS	Acute toxicity - Oral category 4 Skin corrosion/irritation - Skin irritation category 2 Serious eye damage/ Eye irritation category 2A	
Signal Word	WARNING	
Hazards statements (H)	H302 Harmful if swallowed. H315 Causes skin irritation. H319 Causes serious eye irritation.	
Precautionary statements (P)	P264 Wash the areas of the body that have been in contact with the product after handling. P270 Do not eat, drink or smoke when using this product. P280 Wear protective gloves/protective clothing/eye protection/face protection. P302 + P352 IF ON SKIN: Wash with plenty of soap and water. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P321 Specific treatment (see section 4 of the SDS and on this label). P330 Rinse mouth. P332 + P313 If skin irritation occurs: Get medical advice/attention. P337 + P313 If eye irritation persists: Get medical advice/attention. P362 + P364 Take off contaminated clothing and wash it before reuse. P501 Dispose of contents/container in accordance with local / regional / national / international regulations or contact a specialist waste disposal company. P301 + P312 IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.	
PICTOGRAMS		
Other dangers	NFPA (Risk: 0=No risk; 1=Slight; 2=Moderate; 3=Signifiant; 4=Extreme)	
	Health 2 Fire 0 Reactivity 0 Special danger	

SECTION 03 - COMPOSITION/INFORMATION ON INGREDIENTS

Ingrédients (Dénomination chimique / synonymes)	Numéro CAS et tout identificateur unique	Concentration (%)
CHLORURE DE LITHIUM	7447-41-8	<=100%

SECTION 05 - FIREFIGHTING MEASURES

Flammability	No
Ignition conditions	Not flammable or combustible.
Suitable extinguishing media	Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.
Unsuitable extinguishing media	Data not available.
Hazardous combustion / decomposition products	When heated to decomposition, lithium chloride releases toxic vapors of chloride hydrogen.
Special fire and explosion hazards	Lithium chloride reacts violently with bromine trifluoride. May react violently with incompatible products (Ref Section 10).
Special protective equipment and precautions for firefighters	Discard incompatible substances if this can be done without risk. Firefighters should be equipped with standard protective equipment, fireproof clothing, face mask, gloves, protective boots and, where appropriate, self-contained breathing apparatus.

SECTION 06 - ACCIDENTAL RELEASE MEASURES

Methods and materials for containment and cleaning up / Personnel precautions, protective equipment	Evacuate personnel to safe areas. Ensure adequate ventilation. Use personal protective equipment. Pick up with a shovel or broom, taking care not to scatter dust. Clean and rinse with water. Dispose of residues in a container provided for the disposal of hazardous materials. Do not let product enter drains.
--	--

SECTION 07 - HANDLING AND STORAGE

Conditions for safe storage	Hygroscopic. Store in a cool, dry place. Store in a well-ventilated area. Handle and store under inert gas. Keep container tightly closed and store away from heat, moisture, and incompatible products.
Methods of handling	Always open containers slowly to allow any excess pressure to vent. Avoid formation of dust and aerosols. Provide appropriate exhaust ventilation at places where dust or vapor is formed.

SECTION 09 - PHYSICAL AND CHEMICAL PROPERTIES

Physical state	Solid.
Appearance	Poudre cristalline, cristaux ou petites billes de couleur blanche.-
Odour	inodore.
Odour threshold	Data not available
pH	Solution à 5% = pH 5.7 - 7.0 : ce produit est neutre lorsque pur..
Melting point / Freezing point	605°C-
Initial boiling point	1360°C-
Boiling range	Data not available
Flash point	Data not available
Evaporation rate	Data not available
Flammability	No
Lower flammable / Explosive limit	Data not available
Upper flammable / Explosive limit	Data not available
Vapour pressure	1.33 hPa @ 547 °C-
Vapour density	Data not available
Relative density	2.070g/cm ³
Solubility	Très soluble dans l'eau; 569 g/L @ 20 °C. Soluble dans l'alcool et l'acétone.
Partition coefficient water/n-octanol	log Pow: -2.7-
Auto-ignition temperature	Data not available
Decomposition temperature	Data not available
Viscosity	Data not available

SECTION 10 - STABILITY AND REACTIVITY

Reactivity	Non-reactive under normal conditions.
Chemical stability	Stable under recommended storage conditions.
Possibility of hazardous reactions	Stable under normal conditions.
Conditions of instability (Including sensitivity to shock / static discharge / vibration)	Extremely hygroscopic. Absorbs moisture from the air and becomes liquid. Exposure to water. Exposure to light. Incompatible products.
Incompatible material	Strong oxidizing agents (nitric acid, perchloric acid, peroxides, chlorates and perchlorates), strong acids, bromine trifluoride, heat and moisture.
Hazardous decomposition products	Hazardous decomposition products formed under fire conditions. Hydrogen chloride gas, lithium oxides.

SECTION 11 - TOXICOLOGICAL INFORMATION

LITHIUM CHLORIDE

Routes of exposure	Ingestion, inhalation, skin and eyes.
Acute exposition effects / symptoms:	By exposure route below.
- Eyes	Irritation and may cause inflammation of the conjunctiva.
- Skin	Irritation.
- Inhalation	Irritation of the mucous membranes and respiratory tract. Nervous disorders, cough, dyspnea, headache, dizziness, seizures, nausea and vomiting.
Acute toxicity (Ingestion)	Irritation of the mucous membranes. Gastrointestinal disorders, kidney damage, cramps, diarrhea, headache, dizziness, confusion, prostration, tremors, blurred vision, seizures, nausea and vomiting.
Chronic exposure effects / symptoms	Burning sensation, conjunctivitis, kidney damage, nerve disorders, cough, dyspnea, headache, dizziness, confusion, irritability, sweating, salivation, tremors, fatigue, weight loss and loss of appetite, nausea and vomiting. Prolonged exposure may cause reproductive abnormalities in humans.
DL50 (specify species and route of entry)	LD50 Oral - Rat - 526 mg/kg. LD50 Dermal - Rat - > 2000 mg / kg
CL50 (specify species and route of entry)	LC50 Inhalation - Rat - > 5.57 mg / L/ 4h

SECTION 12 - ECOLOGICAL INFORMATION

Ecotoxicity	Toxicity to fish: EC50: 158 mg/L/96h (rainbow trout)
Persistence and degradability	Soluble in water Persistence is unlikely based on information provided.
Bioaccumulative potential	Data not available.
Mobility in soil	Soluble in water. Probable mobility in the environment due to its solubility in water.
Other adverse effects	Do not throw residues in the sewer.

SECTION 13 - DISPOSAL CONSIDERATIONS

Waste Disposal Method	Dispose of contents / container in accordance with local / regional / national / international regulations / or contact a specialist waste disposal company.
Contaminated Packaging	Dispose of as unused product.

SECTION 14 - TRANSPORT INFORMATION

UN Number	N/R
UN Proper shipping name	
Transport hazard class(es)	
Packing group	
Limited quantity index	
ERAP Index	
Special precautions	

SECTION 15 - REGULATORY INFORMATION

WHIMS CANADA	Acute toxicity - Oral category 4 Skin corrosion/irritation - Skin irritation category 2 Serious eye damage/ Eye irritation category 2A
---------------------	--

SECTION 16 - OTHER INFORMATION

Further information

The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. It does not represent any guarantee of the properties of the product. Laboratoire MAT Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.

Last Update: 7/5/2019