



SAFETY DATA SHEET

SECTION 01 - IDENTIFICATION

Product Identifier	PEPSINE (1:3000)
Other identification (Chemical name, Commercial name, Synonymous)	Pepsin A; PEPSINE POWDER
Product code	PP-0411
Chemical formula	-
Molar weight	
Recommended use and Restrictions on use	For laboratory, school, commercial or industrial use. Not for medical or household use. Do not use for medical, food or household purposes.
Supplier	LABORATOIRE MAT 610, rue Adanac Québec Québec G1C 7B7 418-660-8666 Mon-Fri 8h-16h www.labmat.com labmat@labmat.com
Emergency phone	418-660-8666 Mon-Fri 8h-16h CENTRE ANTI-POISON DU QUÉBEC 800-463-5060
Date SDS	2026-06-05

SECTION 02 - HAZARDS IDENTIFICATION

WHIMS CANADA

- Skin corrosion/irritation - Skin irritation - category 2
- Serious eye damage/eye irritation - Eye irritation - category 2A
- Respiratory or skin sensitization - Skin sensitize - category 1
- Respiratory or skin sensitization - Respiratory sensitizer - category 1
- Specific target organ toxicity - Single exposure - category 3

PICTOGRAMS



Signal Word

DANGER

Hazards statements (H)

- Causes skin irritation
- Causes serious eye irritation
- May cause an allergic skin reaction
- May cause allergy or asthma symptoms or breathing difficulties if inhaled
- May cause respiratory irritation

Precautionary statements (P)

- Wash thoroughly after handling.
- IF ON SKIN: Wash with plenty of water
- Specific treatment (see section 4 on the SDS on this label).
- If skin irritation occurs: Get medical attention.
- Take off contaminated clothing and wash before re-use.
- Wear protective gloves (nitrile, butyl, neoprene), protective clothing and eye and face protection.
- IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- If eye irritation persists get medical attention.
- Avoid breathing mists, gases, vapors and other fumes, or the product itself.
- Contaminated work clothing should not be allowed out of the workplace.
- If skin irritation or a rash occurs: Get medical attention.
- Dispose of contents and container in accordance with local, regional and national regulations, or contact a specialist waste disposal company.
- Wear respiratory protection.
- IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing.
- Use only outdoors or in a well-ventilated area.
- Call a POISON CENTER or a doctor or if you feel unwell.
- Store in a well ventilated place. Keep container tightly closed.
- Store locked up.

Other dangers

NFPA (Risk: 0=No risk; 1=Slight; 2=Moderate; 3=Significant; 4=Extreme)

Health 2
Fire 0
Reactivity 0
Special danger

SECTION 03 - COMPOSITION/INFORMATION ON INGREDIENTS

Component	No. CAS	% Weight
Pepsine 1:3000	9001-75-6	95-100%

SECTION 04 - FIRST AID MEASURE

Eye contact	Wash eyes with large amounts of water for at least 15 minutes while holding eyelids apart to rinse eyes. If irritation persists, seek medical attention.
Skin contact	Wash skin with plenty of water for at least 15 minutes. Remove contaminated clothing. If irritation persists, seek medical attention.
Inhalation	Move the unwell person to the fresh air. If breathing is difficult, give oxygen. Consult a physician.
Ingestion	If the person is conscious, give water to drink. Never give anything by mouth to an unconscious person. Consult a physician.
Most important symptoms and effects (acute and delayed)	Main symptoms of high exposure: Skin, eye and respiratory system irritation. Skin sensitizer. Respiratory sensitizer. Allergic reaction. Ref. section 11.
Immediate medical attention and special treatment, if necessary	Treat according to symptoms. Show this sheet to the attending physician.

SECTION 05 - FIREFIGHTING MEASURES

Suitable extinguishing media	Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.
Unsuitable extinguishing media	No limitations on extinguishing agents.
Combustion products	To our knowledge, the products of combustion have not been fully studied.
Specific hazards of the dangerous product	May form combustible dust concentrations in the air. May react violently with incompatible products (Ref Section 10).
Special protective equipment and precautions for firefighters	Discard incompatible substances if this can be done without risk. Firefighters should be equipped with standard protective equipment, fireproof clothing, face mask, gloves, protective boots and, where appropriate, self-contained breathing apparatus.

SECTION 06 - ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency measures	Evacuate personnel to safe areas. When handling, wear appropriate safety equipment (reference Section 8 for protective equipment to be used). Ensure a good ventilation. Use NIOSH cartridge respiratory protection if necessary or for larger spills. Cut off all sources of ignition. Avoid dust formation. Avoid breathing dust.
Methods and materials for containment and cleaning up	Pick up with a shovel or broom, taking care not to scatter dust. Dispose of residues in a container for disposal of hazardous materials. When handling, wear suitable safety equipment. Discharge into the environment must be avoided.

SECTION 07 - HANDLING AND STORAGE

Conditions for safe storage	Store in a cool, dry, and well-ventilated place. Keep container tightly closed and store away from heat, water, moisture, and incompatible products (ref. section 10). Recommended storage temperature: 2-8 ° C.
Methods of handling	Avoid contact with eyes, skin and clothing. Avoid ingestion and inhalation. Wear personal protective equipment (ref. section 8) when handling. Always ensure good ventilation. Apply the usual standard hygiene rules: Wash your hands after use. Do not eat or drink during use.

SECTION 08 - EXPOSURE CONTROLS/PERSONAL PROTECTION

WORKPLACE CONTROL PARAMETERS

Component	CAS-No.	Value
Pepsine 1:3000	9001-75-6	No occupational exposure limits established by region-specific regulators - Quebec, Alberta, Ontario, British Columbia.

Respiratory	If permitted levels are exceeded, use NIOSH cartridge respiratory protection, or an air-supplied respirator.
Gloves	Handle with protective gloves. Suggested material: Nitrile. Butyl. Neoprene. The type, thickness and length of the glove must be chosen according to the use, the concentration of the product, as well as the duration of use. Replace gloves regularly for better protection.
Eyes	Safety goggles with safety shutters.
Shoes	Use safety shoes.
Clothes	Labcoat. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.
Engineering control	Use fan. Have safety showers and eyewash stations in the workplace in case of an emergency and a ventilation system to maintain the level of concentrations in the air below the exposure limit values.

SECTION 09 - PHYSICAL AND CHEMICAL PROPERTIES

Physical state	Solid
Color	White to pale yellow
Odour	Data not available
Odour threshold	Data not available
Melting point and freezing point	Data not available
Boiling point	Data not available
Flammability	Combustible
Lower flammable / Explosive limit	Data not available
Upper flammable / Explosive limit	Data not available
Flash point	Data not available
Auto-ignition temperature	Data not available
Decomposition temperature	Data not available
pH	Data not available
Kinematic viscosity	Not applicable
Solubility	Data not available
Partition coefficient water/n-octanol	Data not available
Vapour pressure	Data not available
Relative density	Data not available
Vapour density	Data not available
Particle characteristics	Powder or crystals

SECTION 10 - STABILITY AND REACTIVITY

Reactivity	May react violently with incompatible substances. Fine dusts of organic materials have the potential to form a combustible or explosive mixture with air if confined in a small space and subjected to an ignition source.
Chemical stability	Stable under recommended storage conditions.
Possibility of hazardous reactions	May react violently with incompatible substances.
Conditions to avoid	Avoid moisture. Avoid contact with incompatible materials and extreme temperatures.
Incompatible materials	Alcohols. Alkalis. Metal salts. Tannins.
Hazardous decomposition products	Carbon oxides. To our knowledge, the products of decomposition have not been fully studied.

SECTION 11 - TOXICOLOGICAL INFORMATION

PEPSINE (1:3000)

Routes of exposure	Ingestion, inhalation, skin and eye contact.
Acute exposition effects / symptoms:	By exposure route below.
- Eyes	Irritation. To our knowledge, the chemical, physical and toxicological properties have not been fully investigated.
- Skin	Irritation. To our knowledge, the chemical, physical and toxicological properties have not been fully investigated.
- Inhalation	May irritate the respiratory system. To our knowledge, the chemical, physical and toxicological properties have not been fully investigated.
Acute toxicity (Ingestion)	Irritation of the esophagus and digestive tract. To our knowledge, the chemical, physical and toxicological properties have not been fully investigated.
Chronic exposure effects / symptoms	Respiratory sensitizer. Skin sensitization. Allergic reactions. To our knowledge, the chemical, physical and toxicological properties have not been fully investigated.
DL50 (specify species and route of entry)	LD50 Oral - Data not available. LD50 Dermal - Data not available.
CL50 (specify species and route of entry)	LC50 Inhalation - Data not available.

SECTION 12 - ECOLOGICAL INFORMATION

PEPSINE (1:3000)

Ecotoxicity	Data not available.
Persistence and degradability	Data not available.
Bioaccumulative potential	Data not available.
Mobility in soil	Data not available.
Other adverse effects	Data not available.

SECTION 13 - DISPOSAL CONSIDERATIONS

Waste Disposal Method	Dispose of contents and container in accordance with local, regional and national regulations, or contact a specialist waste disposal company.
Contaminated Packaging	Dispose of as unused product.

SECTION 14 - TRANSPORT INFORMATION

UN Number	N/R
UN Proper shipping name	
Transport hazard class(es)	
Packing group	
Limited quantity index	
ERAP Index	
Special precautions	

SECTION 15 - REGULATORY INFORMATION

WHIMIS CANADA

- Skin corrosion/irritation - Skin irritation - category 2
- Serious eye damage/eye irritation - Eye irritation - category 2A
- Respiratory or skin sensitization - Skin sensitize - category 1
- Respiratory or skin sensitization - Respiratory sensitizer - category 1
- Specific target organ toxicity - Single exposure - category 3

SECTION 16 - OTHER INFORMATION

CNESST: Commission des normes, de l'équité et de la santé et sécurité au travail

NIH: National institute of health (U.S. National Library of Medicine)

ECHA: Agence Européenne de Chimie

ACGIH : American Conference of Governmental Industrial Hygienists

AIHA : American Industrial Hygiene Association

VECD: Valeur d'exposition courte durée

VEMP: Valeur d'exposition moyenne pondérée

CMRG : Chemical Manufacturer's Recommended Guidelines

OSHA : United States Department of Labor - Occupational Safety and Health Administration

TLV : Threshold limit value

TWA: Time-Weighted-Average

STEL: Short Term Exposure Limit

CEIL: Ceiling

RSST: Règlement sur la santé et sécurité au travail (Québec)

INRS: l'Institut national de recherche et de sécurité pour la prévention des accidents du travail et des maladies professionnelles (France)

The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. It does not represent any guarantee of the properties of the product. Laboratoire MAT Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.

LAST UPDATE: 2026-06-05